

DOCUMENT NAME

Document Title	User Interface Design Document		
Author(s)	Dung, Mai Van		
Role	Team Member		
Date	10/12/2025	File name:	C1SE.40_User_Interface_Design_ver2.0.docx

REVISION HISTORY

Version	Date	Comments	Author
1.0	29/09/2025	Initial Release	All Members
1.1	16/10/2025	Updated	All Members
1.2	15/11/2025	Updated	All Members
1.3	20/11/2025	Updated	All Members
1.4	05/12/2025	Updated	All Members
2.0	10/12/2025	Official	All Members

TABLE OF CONTENTS

1. Introduction.....	1
1.1. Purpose	1
1.2. Scope	1
2.Interface Design.....	1
2.1.Registration & Login.....	1
2.2. Home Page	4
2.3. Weight & BMI Tracking.....	6
2.4. Camera Management	9
2.5. LifeStyle	11
2.7. Settings	14
2.8. Alerts	15
2.9. AI Doctor.....	17
2.10. Chat Screen	19
2.11. Health Tracking.....	20
2.12. Blood Sugar.....	22
2.13. Blood Pressure	24
2.14. Login / Dashboard.....	26
2.15. ElderCare Admin – Dashboard	27
2.16. Analysis Report.....	30
2.17. UI Specification	32
2.18. Approve device registration	34
2.19. Device details	36
2.20. Device management	37
2.21. User details.....	39
2.22. User management.....	40
2.23. Dashboard.....	42

1. Introduction

1.1. Purpose

This specification covers following:

- Documentation assists in achieving a comprehensive understanding of the HealthApp mobile application's user interface and its detailed description and functionality.
- With this documentation, the development team can understand the structural and compositional details to create a consistent and complete health management system.

1.2. Scope

- Overview of the stages and modules of the HealthApp software, including authentication, data tracking, device management, and lifestyle features.
- Create a complete mobile application interface and describe its elements and interactions clearly.

2. Interface Design

2.1. Registration & Login



15:57

HealthApp
Chăm sóc sức khỏe

Đăng ký

Họ và tên

Email

Số điện thoại

Mật khẩu

Xác nhận mật khẩu

Địa chỉ vật lý thiết bị Camera (Tùy ch...

Địa chỉ vật lý thiết bị Vòng tay (Tùy ch...

Đăng ký

Đã có tài khoản? **Đăng nhập**

Registration				
No.	Field name	Type	Require	Target
1	Background	Background	Background color: #101522	Defines the overall visual background of the registration screen.
2	Registration Form	Screen	Heading: “Registration”	Screen that allows users to create a new account.
3	Full Name	Input Text	Placeholder: “Full Name”	Allows the user to enter their full name.
4	Email	Input Text	Placeholder: “Email”	Allows the user to enter an email address for account registration.
5	Phone Number	Input Text	Placeholder: “Phone Number”	Allows the user to enter a phone number.
6	Password	Input Password	Placeholder: “Password”; Eye icon to show/hide password	Allows the user to set an account password securely.
7	Confirm Password	Input Password	Placeholder: “Confirm Password”; Eye icon to show/hide password	Allows the user to re-enter and confirm the password.
8	Camera Device MAC Address	Input Text (Optional)	Placeholder: “Camera device MAC Address (Optional)”	Allows the user to enter the camera device MAC address (optional).
9	Bracelet Device MAC Address	Input Text (Optional)	Placeholder: “Bracelet device MAC Address (Optional)”	Allows the user to enter the bracelet device MAC address (optional).
10	Registration Button	Button	Background color: Green	Submits the registration information to create a new account.
11	Already have an account?	Text Link	Text color: Green	Redirects the user to the Login screen.

15:57



HealthApp
Chăm sóc sức khỏe

Đăng nhập




Đăng nhập

[Quên mật khẩu?](#)

Chưa có tài khoản? [Đăng ký](#)

Login				
No.	Field Name	Type	Require	Target
1	Login Form	Screen	Heading: "Login"	Screen that allows users to access an existing account.
2	Email / Phone Number	Input Text	Placeholder: "Email / Phone Number"	Allows the user to enter an email address or phone number.
3	Password	Input Password	Placeholder: "Password"; Slashed eye icon to hide/show password	Allows the user to enter their account password securely.
4	Login Button	Button	Background color: Green	Authenticates the user and grants access to the application.
5	Forgot Password?	Text Link	Text color: Green	Redirects the user to the password recovery process.
6	Don't have an account?	Text Link	Text color: Green	Redirects the user to the Registration screen.


2.2. Home Page



Home Page				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Home Page”	Identifies the current screen displayed to the user.
2	Current Weather	Text / Icon	Temperature: 26°C; Weather icon: Cloud / Sun / relevant weather icon	Displays current weather information to provide environmental context.
3	Devices Section	Section Heading	Title: “Devices”	Contains status information for all connected smart devices.
4	Camera Device Status	Card	Device name: “Camera”; Date/Time: 29/10/2025, 14:47:26	Displays the latest update time and current operational status of the camera device.
5	Offline / Expired	Alert Badge	Background color: Red;	Alerts the user to critical device

	Alert		Text: “Device is offline”, “Expired”	connection or subscription issues.
6	View Help Link	Text Link	Text: “View help”	Redirects the user to troubleshooting instructions or device support information.
7	Health Diary Section	Section Heading	Title: “Health Diary”	Provides an overview of key health metrics with access to detailed tracking features.
8	Blood Pressure	Metric Card	Value: 100/75 mmHg; Icon: Heart / Pulse wave	Displays blood pressure data and links to detailed blood pressure tracking.
9	Blood Sugar	Metric Card	Value: -- mmol/L; Icon: Finger pricking	Displays blood sugar status and links to detailed blood sugar tracking.
10	Weight & BMI Index	Metric Card	Title: “Weight & BMI Index”; Icon: Weight scale	Provides access to detailed weight and BMI tracking information.
11	Water Reminder	Metric Card	Title: “Water Reminder”; Icon: Water glass	Links to the water intake reminder and tracking feature.
12	AI Doctor Access	Floating Action Button (FAB)	Position: Center of bottom navigation bar; Icon: Exclamation mark in a circle (!)	Provides quick access to the AI Doctor / consultation feature (intended feature).
13	Bottom Navigation Bar	Navigation	Tabs: Home, Tracking, FAB (AI Doctor), Lifestyle, Settings	Serves as the main navigation mechanism of the application.

2.3. Weight & BMI Tracking



Weight and BMI				
No.	Field Name	Type	Require	Target
1	Screen	Screen	Design and requirement details for Weight & BMI Tracking	Displays the interface for monitoring weight and BMI records over time.
2	BMI Chart	Bar Chart	Max: 33.10; Min: 23.70; Average: 25.65	Visualizes historical BMI trends across time periods (T1 to T5).
3	Weight Chart	Bar Chart	Title: “Weight”; Unit: kg	Visualizes historical weight trends over time.
4	Add Record Button	Button	Text: “+ Add Record”; Background color: Light Green	Redirects the user to the weight and height input screen to add a new record.



History Record				
No.	Field Name	Type	Require	Target
1	Total Records	Card	Text: “Total Records”; Value: 10	Displays the total number of recorded weight and BMI entries.
2	Average BMI	Card	Text: “Average BMI”	Displays the average BMI calculated from all recorded entries.
3	Latest Record Overview	Summary Card	Date & Time: 25/10/2025 • 10:16; Weight: 80.0 kg; BMI: 26.1; Status: Overweight	Provides a quick overview of the most recent weight and BMI record.
4	Filter Tabs	Tabs	Options: All, Underweight, Normal, Overweight, Obese	Allows users to filter history records based on BMI classification.
5	Group Header	Text	October 2025	Groups history records by month and year for better readability.
6	History List	List View (per item)	Weight: 80.0 kg / 90.0 kg (blue column); Status: Overweight / Obese (red/orange text); Details:	Displays detailed historical weight and BMI records in a chronological list


			BMI, Height (e.g., 175 cm / 165 cm); Timestamp: 25/10/2025 • Time	format.
--	--	--	---	---------



Add record				
No.	Field Name	Type	Require	Target
1	New Record Screen	Screen	Title: "New Record"	Provides an interface that allows users to manually enter new weight and height data.
2	Weight Input	Input Field	Current value: 60.0 kg; Pencil icon for editing	Allows the user to enter or edit the current body weight.
3	Height Input	Input Field	Current value: 165.0 cm	Allows the user to enter or edit the current body height.
4	BMI Status	Text / Indicator	Status: Normal; BMI value: 22.0	Displays the calculated BMI and its corresponding health classification.
5	BMI Color Scale	Visual Scale	Color-coded BMI scale with indicator arrow pointing to 22.0 (Normal)	Provides a visual reference to help users understand their BMI status.
6	Personal Info	Text	Gender: Male; Age: N/A	Displays basic personal information relevant to

				BMI calculation.
7	Date & Time Header	Header	Header text for date and time selection	Identifies the section used to select the timestamp of the record.
8	Date & Time Picker	Date / Time Picker	Selectable dates (e.g., Tue 21/1, Wed 22, Thu 23) and time values (e.g., 19:33, 20:34)	Allows the user to select the date and time associated with the record.
9	Add Note Button	Button	Text: “Add Note +”	Opens an input field allowing the user to add notes to the record.
10	Save Button	Button	Text: “Save”; Background color: Green	Saves and finalizes the newly entered weight and height record.


2.4. Camera Management



Camera				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Camera”; Subtitle: “Not connected”	Identifies the current camera screen and displays the connection status.
2	Settings Icon	Icon Button	Icon: Gear / Settings	Provides access to camera-specific configuration settings.
3	Video Stream Area	Content Area	Background: Black; Placeholder text: “Camera not	Displays the live video stream when the camera is connected.

			connected”	
4	Connection Guide Text	Text	Instruction: “Add HLS stream link to watch live home camera”	Guides the user on how to connect the camera via streaming link.
5	Guide Button	Button	Text: “Connection Guide”; Background color: Blue	Navigates to a detailed step-by-step camera connection guide.
6	Demo Link Button	Button	Text: “Use demo link”; Background color: White	Loads a demo video stream to allow users to test the camera feature.
7	Connection Status	Text	Text: “Status: Not connected”	Confirms and displays the current camera connection status.
8	Feature Buttons (Top Row)	Icon Buttons	Icons: Expired, Review, Alarm Notification, Multi-channel Screen	Provides quick access to primary camera-related features.
9	Action Icons (Mid-Screen)	Icon Buttons	Icons: Grid view, Refresh	Allows switching layouts or refreshing the camera stream/status.
10	Timeline / Timestamp	Text	Date & Time: 2025-10-13 16:00	Indicates the timestamp of the video segment being viewed (review mode).
11	Control Buttons (Bottom Row)	Icon Buttons	Icons: Guide, Connect, Fullscreen	Provides core camera control actions for the user.

2.5. Lifestyle



Lifestyle				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Lifestyle”	Identifies the content section related to overall well-being and daily lifestyle activities.
2	Content Filter Tabs	Tabs	Options: All, Sleep Better, Health Check	Filters the displayed lifestyle content based on selected category.
3	Sleep Better Section	Video / Article List Section	Title: “Sleep Better”	Provides videos and articles focused on improving sleep quality.
4	Health Check Section	Tool Card Section	Title: “Health Check”	Contains interactive tools that allow users to self-assess their health conditions.
5	Health Tips Section	Card Section	Title: “Health Tips”	Displays short and practical tips to help users maintain a healthy lifestyle.

6	Articles Section	Article List Section	Title: “Articles”	Provides in-depth educational articles on various health-related topics.
7	Cooking Recipes Section	Card Section	Title: “Cooking Recipes”	Offers healthy cooking ideas and nutrition-friendly recipes.


2.6. Water Reminder



Water Reminder				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Water Reminder”	Identifies the screen dedicated to water intake tracking.
2	Current Intake	Large Text	Value: ml (updates dynamically as the user logs water intake)	Displays the total amount of water consumed during the current day.
3	Daily Goal	Text	Title: “Daily Goal:”; Value: 2000 ml; Icon: Water cup	Displays the target daily water intake the user aims to achieve.
4	Last Drink	Text	Title: “Last drink:”;	Displays the most


			Value: 0 ml (updated with the most recent entry)	recent recorded water intake amount.
5	Number of Drinks	Text	Title: “Number of drinks:”; Value: 0 cups; Icon: Water cup	Displays the total number of water intake entries recorded for the day.
6	History & Stats Link	Text Link	Text: “History and Statistics >”	Navigates the user to a detailed history and statistical view of water intake.
7	Progress Indicator	Progress Bar / Text	Value: Dynamic (based on current intake vs. daily goal)	Visually represents the progress toward achieving the daily water intake goal.
8	Decrease Button	Button	Icon: Minus (-)	Decreases the current water intake amount (if applicable).
9	Increase Button	Button	Icon: Plus (+)	Adds a default amount of water or opens an interface to select an intake amount (e.g., 200 ml).
10	Quick Add Button	Button	Value: 200 ml; Icon: Water bottle	Allows the user to quickly log a fixed amount of water (200 ml).

2.7. Settings



Settings				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Settings”; Subtitle: “Manage account and application”	Identifies the screen used for managing user preferences and application settings.
2	Overview Section	Section Heading	Title: “OVERVIEW”	Groups general application-level settings.
3	Notifications	Toggle Switch	Title: “Notifications”; Description: “Receive health notifications and reminders”; Icon: Bell; Toggle: On / Off	Allows the user to enable or disable all application notifications.
4	Language	List Item	Title: “Language”; Current value: “Vietnamese”; Icon: Globe	Allows the user to change the application language.
5	Data Export	List Item	Title: “Data Export”; Description: “Download your health data”; Icon: Book / Export	Initiates the process of exporting and downloading personal health data.
6	Account Section	Section Heading	Title: “ACCOUNT”	Groups account-related actions and settings.
7	Log Out	List Item / Button	Title: “Log out”; Icon: Door	Logs the user out of the current account session securely.

2.8. Alerts



Alerts				
No.	Field name	Type	Require	Target
1	Page Header	Heading	Title: " Latest Alert "	Identifies the section showing the latest and historical alerts.
2	Latest Alert Card	Card		Displays the detailed information for the most recent alert.
3	Video Placeholder/Icon	Icon/Area	Icon: Play button	The intended area to play a video or show an image capturing the alert event.
4	Alert Type Badge	Badge	Text: " Falling Incident " Background: Red	Highlights the type of incident that occurred.
5	Timestamp	Text	Value: 09:30 AM, October 15, 2024	The time the alert was recorded.
6	Location/Details	Text	Text: " Detected fall in the living room "	Provides details about the event and the location of detection.

7	Contact Medical Button	Button	Text: " Contact medical support" Background: Green	An urgent call-to-action button to contact medical services.
8	First Aid Guide Button	Button	Text: " First aid guide " Background: Yellow	A button providing first aid instructions for this type of incident.
9	History Section	Section Heading	Title: " Alert History "	A segment displaying previously occurred alerts.
10	Historical Alert - Fall	List Item/Card	Text: "Sự cố té ngã" Timestamp: 09:30 AM, October 14, 2024 Status: " PROCESSED "	A record of a processed falling incident alert.
11	Historical Alert - Abnormal	List Item/Card	Text: " Abnormal condition " Timestamp: 07:15 AM, October 12, 2024 Status: " PROCESSED "	A record of a processed abnormal health condition alert.
12	Bottom Navigation	Navigation	Tabs: Home, Tracking, Alerts -FAB, Lifestyle, Settings	The main navigation bar for the application.


2.9. AI Doctor



AI Doctor				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “AI Doctor Support”	Identifies the screen dedicated to AI-powered health consultation and support.
2	Chat Input Card	Card / Banner	Text: “Do you have any questions for the doctor?”; Icon: Doctor / Female avatar	Invites the user to enter a health-related question or concern.
3	Ask Question Button	Button	Text: “Ask question”; Icon: Chat bubble	Initiates a consultation or conversation with the AI Doctor.
4	Specialist Section	Section Heading	Title: “Specialist Doctors”	Groups options for selecting a medical specialist.
5	View All Link	Text Link	Text: “See all >”	Redirects the user to a complete list of

				available medical specialists.
6	Cardiologists Option	Icon / Text Label	Text: “Cardiologists”; Icon: Male avatar with heart	Starts a consultation focused on cardiovascular-related issues.
7	Hematologists Option	Icon / Text Label	Text: “Hematologists”; Icon: Female avatar with blood drop	Starts a consultation focused on blood-related conditions.
8	Endocrinologists Option	Icon / Text Label	Text: “Endocrinologists”; Icon: Male avatar with syringe	Starts a consultation focused on endocrine-related conditions.
9	FAQs Section	Section Heading	Title: “Frequently Asked Questions”; Icon: Lightbulb	Displays a list of suggested or common health-related questions.
10	FAQ Item 1	Question Card	Text: “My blood pressure is a bit high, do I have high blood pressure?”; Icon: Chat bubble	Provides a predefined question related to blood pressure concerns.
11	FAQ Item 2	Question Card	Text: “I feel a bit nauseous...”	Provides a predefined question related to nausea symptoms.
12	FAQ Item 3	Question Card	Text: “What could be the cause of my headache?”	Provides a predefined question related to headache causes.
13	FAQ Item 4	Question Card	Text: “Did I catch a cold again? What are the...”	Provides a predefined question related to common cold symptoms.


2.10. Chat Screen



Chat Screen				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Chat”	Identifies the screen as the direct chat and consultation interface.
2	Back Button	Icon Button	Icon: Back arrow	Navigates the user back to the previous screen (AI Doctor Support screen).
3	AI Doctor Message	Chat Bubble (Received)	Sender: “Emma”; Content: “ <i>Hello, my name is Dr. Emma.....</i> ”; Timestamp: 2025-10-06 20:42:40	Displays the initial greeting and consultation prompt from the AI Doctor.
4	Message Options	Icons / Buttons	Icons: Copy, Gallery (Image)	Allows the user to copy message content or open the image gallery related to the conversation.
5	Input Field	Text Input	Placeholder: “My blood pressure is a bit high do I have high blood pressure?”	Allows the user to type a message or health-related question.
6	Send Button	Button	Icon: Paper plane	Sends the typed message to the AI Doctor for processing.

7	Keyboard Options	Icons / Buttons	Icons: Keyboard switch, Gallery, GIF, Document, Settings, Emoji / Sticker	Provides extended input options such as media, documents, emojis, and keyboard settings.
---	------------------	-----------------	---	--


2.11. Health Tracking



Heath Tracking				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Health Tracking”	Identifies the screen as the health metrics overview dashboard.
2	Sync Button	Icon Button	Text: “Sync”; Icon: Refresh / Sync	Allows the user to manually synchronize the latest health data from connected devices.
3	Summary Header	Text	Text: “Today’s Overview”	Serves as the title for the section displaying key health metrics for the current day.
4	Heart Rate Card	Metric Card	Title: “Heart Rate”; Status: Normal; Value: 76 bpm	Displays the current heart rate along with its health status classification.
5	Steps Card	Metric Card	Title: “Steps”; Goal: 8,000 steps; Value: 6,420 / day	Displays the number of steps achieved compared to the daily goal.
6	BMI Card	Metric	Title: “BMI”;	Displays the current BMI

		Card	Status: Normal; Value: 21.3	value and its classification status.
7	Last Sync Status	Card	Title: “Sync”; Details: “Today, 09:42”	Displays the timestamp of the most recent health data synchronization.
8	Weekly Trend Header	Section Heading	Text: “Weekly Trend”	Identifies the section that visualizes health data trends for the week.
9	Weekly Trend Chart	Bar Chart	Data: Steps (Mon–Sun); Average: 6,653 steps/day	Visualizes the weekly trend of the selected health metric (default: steps).
10	Chart Filter Button	Button	Text: “Steps”	Allows the user to switch the metric displayed in the weekly trend chart (e.g., Steps, Heart Rate).
11	Recent Metrics Header	Section Heading	Text: “Recent Metrics”	Serves as the title for the section displaying additional health data.
12	Water Intake Card	List Item / Card	Title: “Water Intake”; Details: “Today”; Value: 1.2 L	Displays the total water intake recorded for the current day.
13	Calorie Consumption Card	List Item / Card	Title: “Calories Consumed”; Details: “Morning”; Value: 530 kcal	Displays calorie consumption for a specific time period.
14	Sleep Card	List Item / Card	Title: “Sleep”; Value: 6h 45m	Displays sleep duration or sleep quality information.
15	Bottom Navigation Bar	Navigation	Tabs: Home, Tracking, FAB (AI Doctor), Lifestyle, Settings	Provides the main navigation structure of the application.


2.12. Blood Sugar



Blood Sugar				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Blood Sugar”	Identifies the screen dedicated to tracking blood sugar metrics.
2	Current Metric Section	Card	Title: “Current Metric”	Displays the most recent blood sugar measurement result.
3	Current Value	Large Text	Value: 90 mg/dL	Shows the latest recorded blood sugar reading.
4	Measurement Unit	Text / Badge	Unit: mg/dL	Indicates the unit of measurement for blood sugar values.
5	Target Range	Badge	Text: “Target: 80–130 mg/dL”	Displays the user’s target blood sugar range.
6	Status	Badge	Text: “Status: Good”	Evaluates the current blood sugar level relative to the target range.
7	Trend	Badge	Text: “Trend: +48 mg/dL”	Indicates the change compared to the previous measurement cycle.
8	Quick Record Section	Card	Title: “Quick Record”	Allows quick entry of commonly used blood sugar values.
9	Quick Value	Button	Value: 80 mg/dL	Quickly logs a blood sugar

	Button 1			value of 80 mg/dL.
10	Quick Value Button 2	Button	Value: 110 mg/dL	Quickly logs a blood sugar value of 110 mg/dL.
11	Quick Value Button 3	Button	Value: 140 mg/dL	Quickly logs a blood sugar value of 140 mg/dL.
12	Quick Value Button 4	Button	Value: 180 mg/dL	Quickly logs a blood sugar value of 180 mg/dL.
13	Add Button	Button	Text: “+ Add”	Opens the detailed or custom blood sugar data entry interface.
14	7-Day Trend Section	Card	Title: “7-Day Trend”	Displays charts and statistics summarizing the last 7 days.
15	Trend Line Chart	Line Chart	—	Visualizes blood sugar fluctuations over the past 7 days.
16	Average Value	Metric	Text: “Average”; Value: 120.2 mg/dL	Displays the average blood sugar level over 7 days.
17	Lowest Value	Metric	Text: “Lowest”; Value: 87 mg/dL	Displays the lowest blood sugar value recorded over 7 days.
18	Highest Value	Metric	Text: “Highest”; Value: 171 mg/dL	Displays the highest blood sugar value recorded over 7 days.
19	Recent History Header	Section Heading	Title: “Recent History”	Identifies the section displaying the most recent detailed blood sugar records.


2.13. Blood Pressure



Blood Pressure				
No.	Field Name	Type	Require	Target
1	Page Header	Heading	Title: “Blood Pressure Tracking”	Identifies the screen dedicated to monitoring blood pressure metrics.
2	Back Button	Icon Button	Icon: Back arrow	Navigates the user back to the previous screen.
3	Sync Button	Icon Button	Text: “Sync”; Icon: Refresh / Sync	Synchronizes the latest blood pressure data from connected devices.
4	Latest Record Time	Text / Icon	Icon: Pulse; Text: 01/12/2025 • 16:48	Displays the timestamp of the most recent blood pressure record.
5	Add Button	Button	Text: “+ Add”; Background color: Blue	Opens the interface for manual entry of a new blood pressure record.
6	Latest SYS Value	Metric Card	Title: “SYS (Systolic)”; Value: 124 mmHg	Displays the most recent systolic blood pressure value.
7	Latest DIA Value	Metric Card	Title: “DIA (Diastolic)”; Value: 79 mmHg	Displays the most recent diastolic blood pressure value.
8	Latest Pulse Value	Metric Card	Title: “Pulse”; Value: 72 bpm	Displays the latest heart rate (pulse) value.
9	MAP Value	Metric Card	Title: “MAP	Displays the calculated

			(Mean Arterial Pressure)”; Value: 94 mmHg	mean arterial pressure value.
10	PP Value	Metric Card	Title: “PP (Pulse Pressure)”; Value: 45 mmHg	Displays the pulse pressure derived from SYS and DIA values.
11	Average Value	Metric Card	Title: “Average”; Value: 126 / 81	Displays the average blood pressure values over a defined period.
12	7-Day Trend Header	Section Heading	Text: “7-Day Trend (SYS)”	Identifies the chart section showing 7-day systolic blood pressure trends.
13	Trend Range	Text	Min: 113; Max: 140	Displays the minimum and maximum systolic values in the 7-day range.
14	Trend Line Chart	Line Chart	—	Visualizes systolic blood pressure (SYS) changes over the last 7 days.
15	Chart Data Points	List / Labels	Date/Time and values (e.g., 16:48 01/12/2025 • 124/79)	Displays detailed data points used to generate the trend chart.
16	History Header	Section Heading	Title: “Measurement History”	Identifies the section listing detailed blood pressure records.
17	Filter Button	Button	Text: “Filter”; Icon: Filter	Opens filtering and sorting options for the measurement history.
18	Latest History Item	List Item	Value: 124/79; Details: MAP 94 • PP 45 • 72 bpm; Status: Hypertension Grade 1; Time: 16:48 01/12/2025	Displays the most recent detailed blood pressure measurement record.

2.14. Login / Dashboard



Login/Dashboard				
No.	Field name	Type	Require	Target
1	Page Container	Layout / Card	Gradient background (Pink → Orange)	Main container holding the login interface.
2	Illustration Section	Image / Card	3D illustration image	Visual representation to enhance user experience and engagement.
3	Welcome Text	Heading	Text: "Welcome back to"	Greets the user returning to the system.
4	Page Title	Heading	Text: "Dashboard"	Identifies the destination after successful login.
5	Login Form Section	Form	Email + Password fields	Area where users enter their authentication information.
6	Email Input	Input (Text)	Placeholder: "Email"	Allows the user to enter their registered email address.
7	Email Icon	Icon	Envelope icon	Visual cue indicating email input field.
8	Password Input	Input (Password)	Placeholder: "Password"	Allows the user to enter their password securely.
9	Password Icon	Icon	Lock icon	Visual cue indicating password input field.
10	Login Button	Button	Text: "Login"	Submits login credentials for authentication.
11	Button Style	UI Style	Rounded, bordered	Enhances button visibility and usability.

12	Card Shadow	UI Effect	Box-shadow enabled	Creates depth and separates the card from background.
13	Color Theme	UI Style	Pink–Orange gradient	Provides a modern and friendly visual appearance.
14	Typography	UI Style	Bold title, normal inputs	Improves readability and visual hierarchy.
15	Responsive Layout	Layout	Center aligned	Ensures proper display across different screen sizes.


2.15. Eldercare Admin – Dashboard

Eldercare Admin				
No.	Field Name	Type	Require	Target
1	Top Navigation Bar	Header	Logo + Title + Search + Icons	Provides global navigation and quick access actions across the system.
2	System Logo	Image	Text: ElderCare Admin	Identifies the system branding.
3	Page Title	Heading	Text: “DASHBOARD”	Indicates the currently active system page.
4	Search Bar	Input (Search)	Placeholder: “Search for something”	Allows quick searching within the admin system.
5	Notification Icon	Icon	Bell icon	Displays system notifications and alerts.
6	Settings Icon	Icon	Gear icon	Navigates to system configuration settings.
7	User Avatar	Avatar	Initials / Profile image	Represents the currently logged-in admin user.
8	Sidebar	Sidebar	Vertical menu	Enables navigation

	Navigation			between main administrative modules.
9	Home Menu Item	Menu Item	Text: “Trang chủ”	Navigates to the main admin dashboard.
10	Account Menu Item	Menu Item	Text: “Tài khoản”	Opens user account and profile management screen.
11	Device Management Menu	Menu Item	Text: “Quản lý thiết bị”	Manages connected elderly care devices.
12	Registration Approval Menu	Menu Item	Text: “Phê duyệt đăng ký”	Handles approval of new user registration requests.
13	Reports Menu	Menu Item	Text: “Báo cáo”	Displays system usage and health-related reports.
14	Settings Menu	Menu Item	Text: “Cài đặt”	Configures system-level preferences.
15	Logout Button	Action Button	Text: “Đăng xuất”	Logs the admin out of the system securely.
16	Profile Card	Card	Avatar + Name + Email	Displays a summary of the admin user’s profile.
17	Profile Avatar	Image / Avatar	Editable	Displays profile image with option to update.
18	User Full Name	Text	Example: “Mai Văn Dũng”	Displays the administrator’s full name.
19	User Email	Text	Unique email address	Displays the registered email address.
20	Join Date	Info Badge	Date & Time	Indicates when the user joined the system.
21	Last Updated Time	Info Badge	Date & Time	Displays the last profile update timestamp.
22	Basic Information Section	Form Section	Title: “Thông tin cơ bản”	Groups editable personal information fields.
23	Full Name Input	Input (Text)	Required	Allows updating the user’s full name.
24	Email Input	Input (Email)	Required – Unique	Allows updating the email address (must be unique).
25	Date of Birth	Input (Date)	Required	Stores the user’s date

				of birth.
26	Current Address	Input (Text)	Optional	Stores the user's current address.
27	Permanent Address	Input (Text)	Optional	Stores the user's permanent residence address.
28	City	Input (Text)	Required	Stores city information.
29	Postal Code	Input (Text / Number)	Optional	Stores postal or ZIP code.
30	Country	Input (Text)	Required	Stores country of residence.
31	Security Section	Section	Title: “Bảo mật”	Groups account security-related actions.
32	Change Password Button	Button	Text: “Đổi mật khẩu”	Allows the user to change the account password.
33	Reset Button	Button	Text: “Đặt lại”	Resets all form fields to their previous values.
34	Save Button	Button	Text: “Lưu”	Saves updated profile and account information.
35	Card Layout Style	UI Style	Rounded corners + Shadow	Enhances readability and provides a modern UI appearance.
36	Responsive Layout	Layout	Desktop-first	Ensures usability across different screen sizes.

2.16. Analysis Report



Analysis Report				
No.	Field name	Type	Require	Target
1	Page Title	Heading	Text: “Báo cáo & Phân tích”	Identifies the analytics and reporting dashboard.
2	Page Description	Text	Overview description	Summarizes user activity, growth trends, and alert system status.
3	Export Summary CSV Button	Button	Text: “Xuất Summary CSV”	Allows exporting summary analytics data in CSV format.
4	Export Full PDF Button	Button	Text: “Xuất PDF toàn trang”	Exports the entire dashboard as a PDF report.
5	Export Top Users CSV Button	Button	Text: “Xuất Top Users CSV”	Exports top active users data in CSV format.
6	Export Combined CSV Button	Button	Text: “Xuất gộp CSV”	Exports combined analytics data into one CSV file.
7	Total Users Card	Statistic Card	Value: 60	Displays the total number of users in the system.
8	Paid Users Card	Statistic Card	Value: 22	Shows the number of paid users.
9	Total Alerts Card	Statistic Card	Value: 29	Displays the total number of alerts generated.
10	Growth Rate Card	Statistic Card	Value: +25%	Indicates the system growth rate over time.
11	User Growth Chart	Bar Chart	Total vs Paid users	Visualizes user growth comparison over time.
12	Chart Legend	Legend	“Tổng”, “Trả phí”	Explains data categories in the growth chart.

13	Alert Trend Chart	Line Chart	Weekly alert data	Shows the trend of alerts by week.
14	Chart Timeline	Axis	Date-based	Displays time progression for charts.
15	User Distribution Chart	Doughnut Chart	Free vs Paid	Visualizes the distribution of free and paid users.
16	Distribution Legend	Legend	Free / Paid	Explains user segmentation.
17	High Activity Users Section	Table Section	Title: “Người dùng hoạt động cao”	Lists users with the highest activity levels.
18	User Name Column	Table Column	Text	Displays the user’s full name.
19	Activity Score Column	Table Column	Percentage (%)	Indicates user activity level.
20	Alert Count Column	Table Column	Number	Shows how many alerts the user generated.
21	Classification Column	Table Column	Badge	Classifies user activity level.
22	Classification Badge – Excellent	Badge	Text: “Xuất sắc”	Indicates excellent user activity.
23	Classification Badge – Positive	Badge	Text: “Tích cực”	Indicates positive user activity.
24	Classification Badge – Average	Badge	Text: “Trung bình”	Indicates average user activity.
25	Table Row	Table Row	Multiple users	Displays detailed analytics per user.
26	Report Date Footer	Text	Date format	Shows the date the report was generated.
27	Card Layout Style	UI Style	Rounded + Shadow	Improves visual clarity and separation.
28	Color Coding	UI Style	Green / Blue / Red / Yellow	Enhances data interpretation through color meaning.
29	Responsive Layout	Layout	Dashboard grid	Ensures proper display on different screen sizes.
30	Data Visualization Library	Chart Component	Bar / Line / Doughnut	Supports visual analytics representation.

2.17. UI Specification

Yêu cầu • e49e8e6f-d75a-11f0-8132-d493900d6c26
Phạm Văn D — d@example.com

Đóng

Serial thiết bị

Loại thiết bị

Người yêu cầu

Email

Số điện thoại

Ghi chú

Vị trí (tùy chọn)

Duyệt


Từ chối

Hủy

UI Specification				
No.	Field name	Type	Require	Target
1	Request Header	Header	Text: “Yêu cầu – {Request ID}”	Identifies the device registration request with a unique ID.
2	Request ID	Text	UUID format	Unique identifier for tracking the request.
3	Requester Info	Text	Name + Email	Displays requester’s basic identification.
4	Close Button	Button	Text: “Đóng”	Closes the request detail view.
5	Device Serial Number	Input (Read-only)	Example: SN-D-0004	Displays the serial number of the requested device.
6	Device Type	Input (Read-only)	Example: JEWELRY	Indicates the type/category of the device.
7	Requester Name	Input (Read-only)	Example: Phạm Văn D	Shows the name of the user who submitted the request.
8	Email	Input (Read-only)	Example: d@example.com	Displays the requester’s email address.
9	Phone Number	Input (Read-only)	Example: 0901000004	Shows the requester’s contact phone number.
10	Note	Textarea (Read-only)	Text: “Chờ duyệt”	Displays the current status or note of the request.
11	Location (Optional)	Input (Read-only)	Example: Phòng khách	Indicates the intended device location (if provided).
12	Approval Button	Action Button	Text: “Duyệt”	Approves the device registration request.

13	Reject Button	Action Button	Text: “Từ chối”	Rejects the request with administrative decision.
14	Cancel Button	Action Button	Text: “Hủy”	Cancels the current approval action.
15	Button Color – Approve	UI Style	Green	Indicates positive confirmation action.
16	Button Color – Reject	UI Style	Red	Indicates denial or negative action.
17	Form Layout	Layout	Two-column form	Organizes information clearly and efficiently.
18	Read-only Fields	Field State	Disabled	Prevents editing of request data by admin.
19	Status Handling	Logic	Pending / Approved / Rejected	Manages the lifecycle of the request.
20	Validation	System Rule	Required fields present	Ensures request data completeness before action.
21	Confirmation Action	System Action	Popup / Toast	Confirms approval or rejection result.
22	Audit Log	Backend Logic	Admin action recorded	Logs admin decisions for traceability.
23	Security Constraint	Access Control	Admin only	Restricts access to authorized administrators.
24	Responsive Layout	Layout	Desktop-first	Ensures consistent UI on different screen sizes.

2.18. Approve device registration



Approve device registration				
No.	Field name	Type	Require	Target
1	Page Title	Heading	Text: “Phê duyệt đăng ký thiết bị”	Identifies the device registration approval management screen.
2	Page Description	Text	Instruction text	Explains the purpose of approving or rejecting device registration requests.
3	Refresh Button	Button	Text: “Làm mới”	Reloads and refreshes the request list data.
4	Search Bar	Input (Search)	Placeholder: “Tìm theo thiết bị, người dùng, email...”	Allows searching requests by device, user name, or email.
5	Status Filter	Dropdown	Default: “Tất cả”	Filters requests by status (All / Pending / Approved / Rejected).
6	Request Table	Data Table	Multiple rows	Displays all device registration requests.
7	Request ID Column	Table Column	UUID	Displays the unique request identifier.
8	Device Column	Table Column	Serial + Type	Shows device serial number and device type.
9	User Column	Table Column	Full name	Displays the requesting user’s name.

10	Email Column	Table Column	Email	Displays the requester's email address.
11	Created Time Column	Table Column	Date & Time	Shows when the request was created.
12	Status Column	Table Column	Badge	Indicates the current request status.
13	Status Badge – Pending	Badge	Text: “Đang chờ”	Indicates request is waiting for approval.
14	Status Badge – Approved	Badge	Text: “Đã duyệt”	Indicates request has been approved.
15	Action Column	Table Column	Action Icons	Contains available actions for each request.
16	View Detail Action	Icon Button	Eye icon	Opens detailed view of the request.
17	Approve Action	Icon Button	Check icon	Approves the device registration request.
18	Reject Action	Icon Button	Cross icon	Rejects the device registration request.
19	Pagination Info	Text	“Hiển thị 1–3 trên 3”	Displays current result range and total items.
20	Pagination Control	Pagination	Page navigation	Allows navigation between pages of requests.
21	Table Row Highlight	UI Style	Hover effect	Improves readability and user interaction.
22	Badge Color – Pending	UI Style	Yellow	Visual indication for pending status.
23	Badge Color – Approved	UI Style	Green	Visual indication for approved status.
24	Confirmation Handling	System Logic	Modal / Toast	Confirms approve or reject actions.
25	Access Control	Security Rule	Admin only	Restricts access to authorized administrators.
26	Audit Logging	Backend Logic	Action recorded	Records approval/rejection actions for traceability.
27	Responsive Layout	Layout	Table responsive	Ensures usability on different screen sizes.

2.19. Device details

Thiết bị • SN-L-0012

e49de973-d75a-11f0-8132-d493900d8c26

Serial thiết bị

SN-L-0012

Loại thiết bị

JEWELRY

Vị trí

Trạng thái

Ngoại tuyến

Hoạt động gần nhất

Lưu

Hủy

Device details				
No.	Field name	Type	Require	Target
1	Page Header	Header	Text: “Thiết bị – {Serial Number}”	Identifies the device detail screen with its serial number.
2	Device ID	Text	UUID	Unique internal identifier of the device.
3	Close Button	Button	Text: “Đóng”	Closes the device detail view.
4	Device Serial Number	Input (Read-only)	Example: SN-L-0012	Displays the unique serial number of the device.
5	Device Type	Input (Read-only)	Example: JEWELRY	Indicates the category/type of the device.
6	Location	Input (Text)	Optional	Allows entering or updating the device location.
7	Device Status	Dropdown	Values: Online / Offline	Displays and updates the current device status.
8	Status Value – Offline	Dropdown Option	Text: “Ngoại tuyến”	Indicates the device is currently offline.
9	Last Activity	Input (Read-only)	Date & Time	Displays the most recent activity time of the device.
10	Save Button	Primary Button	Text: “Lưu”	Saves any changes made to the device information.
11	Cancel Button	Secondary Button	Text: “Hủy”	Cancels changes and restores previous values.
12	Form Layout	Layout	Two-column layout	Organizes device information clearly.
13	Editable Fields	Field State	Location, Status	Specifies which fields can be edited.
14	Read-only Fields	Field State	Serial, Type, Activity	Prevents modification of critical device data.
15	Validation Rule	System Rule	Required format	Ensures valid data before saving changes.
16	Update Confirmation	System Action	Toast / Notification	Confirms successful update of device information.

17	Error Handling	System Action	Validation message	Displays errors if update fails.
18	Audit Logging	Backend Logic	Update recorded	Logs device updates for traceability.
19	Access Control	Security Rule	Admin only	Restricts access to authorized administrators.
20	Responsive Design	Layout	Desktop-first	Ensures usability across different screen sizes.

2.20. Device management

Quản lý thiết bị
Thiết bị, gắn thiết bị và trạng thái

Tổng thiết bị: 7 Trực tuyến: 4 Đang chờ: 1

Q. Tìm kiếm: serial, id, loại, vị trí... Tất cả 7 kết quả

Serial	Loại	Chủ sở hữu	Vị trí	Hoạt động gần nhất	Trạng thái	Thao tác
SN-L-0012 e4992b80-d75a-11f0-8132-d493900d6c26	JEWELRY	e4992b80-d75a-11f0-8132-d493900d6c26	—	—	Ngoại tuyến	
SN-K-0011 e4990780-d75a-11f0-8132-d493900d6c26	JEWELRY	e4990780-d75a-11f0-8132-d493900d6c26	—	—	Ngoại tuyến	
SN-I-0009 e4982a23-d75a-11f0-8132-d493900d6c26	JEWELRY	e4982a23-d75a-11f0-8132-d493900d6c26	—	—	Đang chờ	
SN-H-0008 e4982061-d75a-11f0-8132-d493900d6c26	JEWELRY	e4982061-d75a-11f0-8132-d493900d6c26	—	—	Trực tuyến	
SN-D-0004 e4977d9a-d75a-11f0-8132-d493900d6c26	JEWELRY	e4977d9a-d75a-11f0-8132-d493900d6c26	Ban công	3/29/2025, 12:25:00 PM	Trực tuyến	
SN-B-0002 e4962676-d75a-11f0-8132-d493900d6c26	JEWELRY	e4962676-d75a-11f0-8132-d493900d6c26	Phòng ngủ	2/15/2025, 10:45:00 AM	Trực tuyến	
SN-A-0001 e4962676-d75a-11f0-8132-d493900d6c26	JEWELRY	e4962676-d75a-11f0-8132-d493900d6c26	Phòng khách	1/11/2025, 9:35:00 AM	Trực tuyến	

Hiển thị 1-7 trên 7 < Trang 1 / 1 >

Device management				
No.	Field name	Type	Require	Target
1	Page Title	Heading	Text: “Quản lý thiết bị”	Identifies the device management screen.
2	Page Description	Text	Instruction text	Describes device list, assignment, and status management.
3	Refresh Button	Button	Text: “Làm mới”	Reloads the device list and statistics.
4	Total Devices Card	Statistic Card	Value: 7	Displays total number of devices in the system.
5	Online Devices Card	Statistic Card	Value: 4	Shows number of devices currently online.
6	Pending Devices Card	Statistic Card	Value: 1	Shows number of devices pending approval or activation.
7	Search Bar	Input (Search)	Placeholder: “Tìm kiếm: serial, id, loại, vị trí...”	Allows searching devices by serial, ID, type, or location.
8	Status Filter	Dropdown	Default: “Tất cả”	Filters devices by status (All / Online / Offline /

				Pending).
9	Result Count	Text	Example: “7 kết quả”	Displays the number of matching devices.
10	Device Table	Data Table	Multiple rows	Displays the list of devices.
11	Serial Column	Table Column	Serial + Device ID	Displays the device serial number and internal ID.
12	Device Type Column	Table Column	Example: JEWELRY	Indicates the type/category of the device.
13	Owner Column	Table Column	User ID	Displays the owner (user) of the device.
14	Location Column	Table Column	Text	Shows the assigned physical location of the device.
15	Last Activity Column	Table Column	Date & Time	Displays the most recent activity timestamp.
16	Status Column	Table Column	Badge	Shows current device status.
17	Status Badge – Online	Badge	Text: “Trực tuyến”	Indicates the device is currently online.
18	Status Badge – Offline	Badge	Text: “Ngoại tuyến”	Indicates the device is offline.
19	Status Badge – Pending	Badge	Text: “Đang chờ”	Indicates the device is pending approval or setup.
20	Action Column	Table Column	Icon buttons	Contains available actions for each device.
21	View Detail Action	Icon Button	Eye icon	Opens device detail view.
22	Delete Device Action	Icon Button	Trash icon	Removes the device from the system.
23	Pagination Info	Text	“Hiển thị 1–7 trên 7”	Displays current result range and total items.
24	Pagination Control	Pagination	Page navigation	Allows navigation between pages if applicable.
25	Table Row Hover	UI Style	Hover effect	Improves readability and interaction.
26	Color Coding	UI Style	Green / Red / Yellow	Enhances visual distinction of device statuses.
27	Confirmation Dialog	System Logic	Modal	Confirms destructive actions such as deletion.
28	Audit Logging	Backend Logic	Action recorded	Logs device changes for traceability.

29	Access Control	Security Rule	Admin only	Restricts access to device management to admins.
30	Responsive Layout	Layout	Dashboard grid	Ensures usability across different screen sizes.

2.21. User details

Người dùng • Demo Growth M0 Free 1
7710557e-d75c-11f0-8132-0493900dc26

Họ và tên: Demo Growth M0 Free 1 Email: demo_growth_m0_f1@example.com


Số điện thoại: 0907000001 Trạng thái: Hoạt động

[Lưu](#) [Hủy](#)

User details				
No.	Field name	Type	Require	Target
1	Page Header	Header	Text: “Người dùng – {User Name}”	Identifies the user detail screen with user name.
2	User ID	Text	UUID	Unique internal identifier of the user.
3	Close Button	Button	Text: “Đóng”	Closes the user detail view.
4	Full Name	Input (Text)	Required	Displays and allows editing of the user’s full name.
5	Email	Input (Email, Read-only)	Required – Unique	Displays the user’s registered email address.
6	Phone Number	Input (Text)	Required	Displays and allows editing of the user’s phone number.
7	User Status	Dropdown	Values: Active / Inactive	Displays and updates the user’s current status.
8	Status Value – Active	Dropdown Option	Text: “Hoạt động”	Indicates the user account is active.
9	Save Button	Primary Button	Text: “Lưu”	Saves changes made to user information.
10	Cancel Button	Secondary Button	Text: “Hủy”	Cancels changes and restores previous values.
11	Form Layout	Layout	Two-column layout	Organizes user information clearly and consistently.

12	Editable Fields	Field State	Name, Phone, Status	Specifies which user attributes can be updated.
13	Read-only Fields	Field State	Email, User ID	Prevents modification of critical user identifiers.
14	Validation Rule	System Rule	Required fields, format check	Ensures valid input before saving changes.
15	Update Confirmation	System Action	Toast / Notification	Confirms successful update of user information.
16	Error Handling	System Action	Validation message	Displays errors if update fails.
17	Audit Logging	Backend Logic	Update recorded	Logs user updates for traceability.
18	Access Control	Security Rule	Admin only	Restricts access to user management functions.
19	Data Privacy	Security Rule	Mask sensitive data	Ensures compliance with data protection policies.
20	Responsive Layout	Layout	Desktop-first	Ensures usability across different screen sizes.

2.22. User management




User management				
No.	Field name	Type	Require	Target
1	Page Title	Heading	Text: “Quản lý thiết bị”	Identifies the device management screen.

2	Page Description	Text	Instruction text	Describes device list, assignment, and status management.
3	Refresh Button	Button	Text: “Làm mới”	Reloads the device list and statistics.
4	Total Devices Card	Statistic Card	Value: 7	Displays total number of devices in the system.
5	Online Devices Card	Statistic Card	Value: 4	Shows number of devices currently online.
6	Pending Devices Card	Statistic Card	Value: 1	Shows number of devices pending approval or activation.
7	Search Bar	Input (Search)	Placeholder: “Tìm kiếm: serial, id, loại, vị trí...”	Allows searching devices by serial, ID, type, or location.
8	Status Filter	Dropdown	Default: “Tất cả”	Filters devices by status (All / Online / Offline / Pending).
9	Result Count	Text	Example: “7 kết quả”	Displays the number of matching devices.
10	Device Table	Data Table	Multiple rows	Displays the list of devices.
11	Serial Column	Table Column	Serial + Device ID	Displays the device serial number and internal ID.
12	Device Type Column	Table Column	Example: JEWELRY	Indicates the type/category of the device.
13	Owner Column	Table Column	User ID	Displays the owner (user) of the device.
14	Location Column	Table Column	Text	Shows the assigned physical location of the device.
15	Last Activity Column	Table Column	Date & Time	Displays the most recent activity timestamp.
16	Status Column	Table Column	Badge	Shows current device status.
17	Status Badge – Online	Badge	Text: “Trực tuyến”	Indicates the device is currently online.
18	Status Badge – Offline	Badge	Text: “Ngoại tuyến”	Indicates the device is offline.
19	Status Badge – Pending	Badge	Text: “Đang chờ”	Indicates the device is pending approval or setup.
20	Action Column	Table Column	Icon buttons	Contains available actions for each device.
21	View Detail Action	Icon Button	Eye icon	Opens device detail view.

22	Delete Device Action	Icon Button	Trash icon	Removes the device from the system.
23	Pagination Info	Text	“Hiển thị 1–7 trên 7”	Displays current result range and total items.
24	Pagination Control	Pagination	Page navigation	Allows navigation between pages if applicable.
25	Table Row Hover	UI Style	Hover effect	Improves readability and interaction.
26	Color Coding	UI Style	Green / Red / Yellow	Enhances visual distinction of device statuses.
27	Confirmation Dialog	System Logic	Modal	Confirms destructive actions such as deletion.
28	Audit Logging	Backend Logic	Action recorded	Logs device changes for traceability.
29	Access Control	Security Rule	Admin only	Restricts access to device management to admins.
30	Responsive Layout	Layout	Dashboard grid	Ensures usability across different screen sizes.

2.23. Dashboard



Dashboard				
No.	Field name	Type	Require	Target
1	Page Title	Heading	Text: “DASHBOARD”	Identifies the main overview dashboard screen.
2	Welcome Banner	Banner / Card	Greeting message	Displays a personalized welcome message to the admin.
3	Welcome Text	Text	“Xin chào {Admin Name}”	Greets the logged-in administrator.
4	Welcome	Text	Inspirational message	Encourages positive

	Description			engagement with the system.
5	Refresh Button	Button	Text: “Làm mới”	Refreshes dashboard data.
6	New Customers Section	Card	Title: “Khách hàng mới”	Displays recently registered users.
7	New Customer Avatar	Avatar	Initials	Represents new users visually.
8	New Customer Status	Badge	Active / Inactive	Shows current status of new users.
9	View All Customers Link	Link	Text: “Xem tất cả”	Navigates to full user list.
10	Quick Action Input	Input	Placeholder: “Enter amount...”	Allows quick data entry or action.
11	Send Button	Button	Text: “Gửi”	Submits quick action input.
12	Total Customers Card	Statistic Card	Value: 60	Displays total number of customers/users.
13	Customer Distribution Chart	Doughnut Chart	New / Active / Inactive	Visualizes customer status distribution.
14	Chart Legend	Legend	Percentage labels	Explains chart segments.
15	Revenue Growth Section	Chart Card	Title: “Tăng Trưởng Doanh Thu”	Displays revenue growth trends.
16	Revenue Bar Chart	Bar Chart	Monthly revenue data	Visualizes revenue changes by month.
17	Revenue Peak Label	Text	Highest revenue note	Highlights peak revenue month and value.
18	Calendar Widget	Calendar	Month view	Displays current month calendar.
19	Calendar Navigation	Button	Previous / Next	Navigates between months.
20	Add Reminder Button	Button	Text: “+ Thêm Lỗi Nhắc”	Adds a new reminder.
21	Upcoming Reminder Card	Card	Reminder details	Displays upcoming reminder information.
22	Reminder Title	Text	Example: “Lỗi nhắc sắp tới”	Identifies reminder section.
23	Reminder Action Link	Link	Text: “Xem chi tiết”	Opens reminder detail view.
24	Alert Statistics Section	Chart Card	Title: “Số lần cảnh báo”	Displays alert frequency statistics.
25	Alert Bar Chart	Bar Chart	Daily alert data	Visualizes number of alerts per day.

26	Alert Legend	Legend	Alert / No Alert	Explains alert chart data.
27	Sidebar Navigation	Sidebar	Menu items	Enables navigation to other modules.
28	Top Navigation Bar	Header	Search + Icons	Provides global search and quick actions.
29	User Profile Icon	Avatar	Initials	Represents logged-in admin.
30	Notification Icon	Icon	Bell	Displays system notifications.
31	Settings Icon	Icon	Gear	Opens system settings.
32	Responsive Layout	Layout	Dashboard grid	Ensures usability across devices.
33	Color Theme	UI Style	Purple / Blue gradient	Enhances visual appeal and hierarchy.
34	Data Refresh Logic	System Logic	Real-time / periodic	Keeps dashboard data up to date.